



GEM Rotator Design Progress


Chandika Annasiwatta
MOLLER collaboration meeting
May 06, 2023

Outline


- Subsystem overview
- Radial motion
- Rotational motion
- Simulations
- Stand adjustors
- Cable arrangements
- Assembly procedure
- Summary

Subsystem Overview

- Rotator to support 28 GEM modules (4 layers of 7 modules) and 2 layers of trigger scintillators (2 layers of 7 scintillators) for Counting Mode measurements
- Must rotate system to cover full azimuth in sequential measurements (Rotation of $360 \text{ degrees} / 7 = 51.4 \text{ degrees}$, with at least three stopping positions)
- Must move radially to extract GEMs and Trigger scintillators fully from scattered electron flux for main Integrating Mode measurements
- Minimize mass in the path of scattered electrons
- Minimize use of ferrous materials




GEM Rotator – Current design




- Three 2.65 m diameter Al wheels
- Rotational motion via chain drive, sprockets on center wheel. Two rollers per wheel.
- V-groove on front wheel, simple roller on back wheel.
- Al (Bosch 80/20 6105-T5) support stand: minimizing Ferrous materials. Bolted not welded.
- Rotational motion: 2-sector (103°) swing, 2 intermediate stops to cover full azimuth.
- Total mass – ~2100 kg
- Stand – ~460 kg

GEM Rotator – Radial motion


Measurement position



Parked position




GEM Rotator – Radial motion




- Radial motion: 14 stepper-motor linear actuators.
- Each moves pair of GEMs and a single trigger scintillator along non-magnetic linear guide rails (Franke FDD guides).
- GEM pair + Scintillator \approx 25 lbs
- Trapezoidal frames supporting the linear motors and guide rails: 6061-T6 aluminum.
- Stepper Motor Linear Actuator : stainless steel screw.
- Radial motion stroke: 760 mm (GEM layers 1&2) and 720 mm (GEM layers 3&4), limited by hard stops at both extremes.
- The 14 linear actuators operate independently:
 - allows insertion of single pair of GEM detectors in a particular semi-septant, for study of any backgrounds in the main detectors caused by the tracking system.
- The holder frame for the GEM and scintillator modules: aluminum & carbon fiber.

GEM Rotator – Rotational motion


- Chain drive
- Max chain tension to hold max unbalance load (GEM holder + 2 GEMs + scintillator)
 - Measurement position \sim 40 lbs
 - Parked position \sim 70 lbs
- ANSI Number 50-SS, 316 Stainless Steel, Working load 150 lbs






GEM Rotator – FEA displacements



Largest displacement:
< 0.7 mm
- Meets specs (< 3 mm)
(was ~1.7 mm in earlier design)





GEM Rotator – FEA Stresses



- 6061-T6 Al - Tensile Strength 3.10e+8 N/m²
- 6105-T5 Al - Tensile Strength 2.21e+8 N/m²
- Alloy steel - Tensile Strength – 7.23e+8 N/m²


GEM Rotator – FEA Stresses



- 6061-T6 Al - Tensile Strength 3.10e+8 N/m²
- 6105-T5 Al - Tensile Strength 2.21e+8 N/m²
- Alloy steel -Tensile Strength 7.23e+8 N/m²

GEM Rotator – Hertzian stress at roller contacts

- Hertzian stress analysis at roller contacts:
 - Flat roller
 - Max contact pressure 141 Mpa
 - V roller
 - Max contact pressure 170 Mpa
- For 6061–T6 Al: fatigue limit: 10^5 repetitions expect < 200 repetitions



0.1 g seismic analysis

- 0.1 g seismic analysis worst cases

- Seismic forces in +x (beam direction): deflection 0.9534 mm in +x (0.9023 mm parked)
- Seismic forces in +x (beam direction): deflection 0.5422 mm in -y
- Seismic forces in -y (downward): deflection 0.5622 mm in y (0.5496 mm parked)


() - Parked	Max Stress on Base leg – (Alloy Steel) (MN/m ²)	Max Stress on Roller support (Al) (MN/m ²)	X (mm)	Y (mm)	Z (mm)
+x (beam direction)	303.8 (306.7)	77.03 (65.56)	0.9534 (0.9023)	-0.5422 (-0.4645)	0.1869 (-0.1883)
-x	326.9 (325.7)	57.44 (51.25)	-0.4482 (-0.5612)	-0.4044 (-0.3489)	0.1472 (-0.2201)
+z	366.3 (360.7)	72.53 (60.28)	0.5115 (0.5029)	-0.5204 (-0.5344)	0.2836 (0.3509)
-z	205.0 (197.6)	65.51 (51.96)	0.5107 (0.4964)	-0.4743 (-0.4272)	-0.2403 (-0.2707)
+y	250.2 (243.0)	60.50 (50.10)	0.4600 (0.4497)	-0.4260 (-0.3660)	0.1498 (-0.1838)
-y (Downward)	305.8 (297.0)	73.95 (61.23)	0.5622 (0.5496)	-0.5207 (-0.4474)	0.1830 (-0.2246)

- 6061-T6 Al - Tensile Strength 310 MN/m²

- Alloy steel -Tensile Strength 723 MN/m²

Buckling analysis



- $1 < \text{BLF}$, Buckling not predicted. The applied loads are less than the estimated critical loads. Buckling is not expected.
- $\text{BLF} < -1$, Buckling not predicted. Buckling is not expected even if you reverse all loads.



Mode No.	Buckling Factor of Safety
1	-1,148.2
2	-1,078.6
3	-534.41
4	-507.13
5	786.67
6	838.95
7	935.26
8	979.26
9	1,015.5
10	1,024.7
11	1,038
12	1,115.9
13	1,137.1
14	1,171.5
15	1,199.6

GEM Rotator – Modal analysis

- Lowest resonant frequency is 12.201 Hz > 10 Hz target value



List Modes



Study name:Modal

Mode No.	Frequency(Rad/sec)	Frequency(Hertz)	Period(Seconds)
1	76.658	12.201	0.081964
2	80.201	12.764	0.078343
3	135.5	21.565	0.046371
4	158.89	25.288	0.039544
5	163.35	25.998	0.038465
6	163.91	26.086	0.038334
7	164.12	26.12	0.038285
8	164.82	26.231	0.038122
9	164.98	26.257	0.038086
10	165.41	26.326	0.037986
11	165.65	26.364	0.037931
12	165.99	26.418	0.037853
13	166.43	26.488	0.037753
14	167.38	26.639	0.037538
15	167.83	26.711	0.037438


Close Save Help

GEM Rotator – Modal analysis


- Lowest resonant frequency is 11.751 Hz > 10 Hz target value



GEM Rotator - Stand adjustors




- Two at the upstream side, and one at the downstream side, grouted to Hall A floor
- ± 15 mm adjustment in the y direction (vertical) and ± 7.8 mm in the x direction. Give Pitch and Yaw adjustment.
- Location of the Rotator in z (beam direction) not critical.
- Roll angle for the GEMs is adjustable by the rotation mechanism.
- The base-leg adjusters follow designs commonly used in Hall A at JLab for adjustable mounting of heavy equipment.






Cable arrangements


- Cables per GEM module
 - 4 HDMI cables (OD 7mm).
 - 1 gas lines (OD 6.35mm).
 - 1 High Voltage(HV) cable (OD 6mm).
- Cables per Trigger scintillator
 - 1 HV cable (OD 6mm).
 - 1 signal cable (OD 6mm).
- Radial motion
 - Option 1: Cable and hose carriers (Loop Ht. 14.4" or 22.2") - Nylon
 - Option 2: Wire bundle attached on two points
- Rotational motion
 - Cable System suspend from above



Assembly procedure



Assembly procedure



Summary

- GEM rotator
 - Mechanical design is mostly done.
 - Analysis (Static, seismic, Modal, Buckling) are completed.
 - Need to add connection details to the rotator wheel and complete the drawing package.
 - Support stand use commercially available parts other than base leg, side plates and bearing support.

- Radial and rotational motion
 - Stepper motor linear actuator is selected.
 - Radial motion motor and chain selected.
 - Need to select the drive controllers.

- Need to finalize the cable arrangements